



ORTHOPAEDIC SUPPORT FOR ACTIVE PEOPLE

- Weak and Unstable Knee joints
- Pain below and or above knee caps



*SUPPORT * ANTI-SLIP
*FLEXIBLE * STABILITY

Adaptive Gel pad to stabilize Patella and Patellar tendon.

Applies just the right Pressure along the Meniscus line.

Added Support along the knee joints.

Aids in Recovery and Pain management.



When my son tore his ACL ten years ago, I had to give him a knee brace to stabilise it. However, as I listened to his comments over the years, I drew similarities with feedback from other patients, describing the changes they would've preferred in current options. By collaborating with a mechanical engineer, we produced a brace design that is tailored for the Asian knee anatomy.

Adapting the brace to the less tapered build of Asian legs allowed us to produce a knee brace which is more form fitting and flexible.

Combined with the elastic material and unique cross-stitch pattern, this brace offers improved tension and stability across the knee, without compromising on comfort

Kneecap instability in Asians also means we suffer more frequent injuries to the area, leading to patella-femoral pain. By reinforcing the soft gel bracing around the kneecap to fix it in place, this brace helps protect against injury and pain.

We also focused on making the brace suitable for the regional climate since the increased temperature and humidity make braces uncomfortable to wear. Our design allowed us to keep a short brace length to keep the area cool, while silicone beads keep the brace from sliding down during intense activity.



Dr Jeffrey Chew Tec Hock
Managing Director

Centre for Robotic Joint Replacement
Consultant Orthopaedic Surgeon
Joint Replacement Specialist & Sports Surgeon

MBBS (Singapore), FRCS (Glasgow), M.Med (Surgery)(Singapore), FICS (USA), FRCS (Edinburgh), FAMS (Orthopaedics), Mayo Alumnus, Mayo Clinic (USA)

Website: www.cfo.com.sg

Email: jeffrey@cfo.com.sg

#10-33/34/35 Mount Elizabeth Novena Specialist Centre, 38 Irrawaddy Road Singapore 329563

Tel: +65 6684 5828 Fax: +65 6684 5829

#05/06/07 Parkway East Medical Centre 319 Joo Chiat Place, Singapore 427989

Tel: +65 6348 2628 Fax: +65 6348 2605

#15-06/07 Connexion Farrer Park Medical Centre 1 Farrer Park Station Road, Singapore 217562

Tel: +65 6738 2628 Fax: +65 6738 2629



Distributed by Delsson Singapore Pte Ltd
Technicalinfo@delsson.com